

 **biocare**



www.biocare.com.cn

iE 101 / iE 300

Digital Electrocardiograph



iE 101

Single Channel Digital
Electrocardiograph




50
mm

iE 300

3-channel Digital
Electrocardiograph



80
mm



Accuracy and reliability in ECG interpretation

- Biocare CardioPro™ ECG analysis program, proven by CSE/AHA/MIT database
- 12-lead simultaneous acquisition
- Complete digital filters and self-adapting filtering
- Up to 300 seconds of rhythm analysis for easier in arrhythmia locating
- Applied the latest Minnesota code

Streamlined workflow, new operating experience

- "One-touch" functional shortcut keys
- Onscreen operation guide to check the ECG procedure
- Multi language voice reminder and visual alarming for lead-off, low battery and no paper
- Support general settings according to different situations for physicians and patients
- 3 report formats: simple report, median beat report and manual report

Data management solutions


- Massive local storage capacity up to 1500 files
- ECG management software (optional)



- About 1.3kg, light in weight



- 5.0-inch high resolution LCD screen



- Built-in rechargeable lithium battery for 4h continuous work



- Peripheral ports: USB, LAN port

Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

Headquarters: 7/F, Block A, Shenzhen Mingyou Industrial Products Exhibition & Procurement Center, Baoyuan Road, Xixiang Sub-district, Bao'an District, 518102 Shenzhen, P.R. China

Factory: 2/F West, 4th Block, Dayang Road South, Fuyong Sub-district, Bao'an District, 518103 Shenzhen, P.R. China

Tel: +86-755-3661 5333

Fax: +86-755-2796 0643

E-mail: sales@biocare.com.cn